

Patient [REDACTED]

Implant Team [REDACTED]

Hospital [REDACTED]

Date 11/24/2022

PATIENT DEMOGRAPHICS (provided by hospital)

| | | | | |
|---|---|--|------------|---------------------|
| Age: 69 | Sex: Male | Height: m | Weight: kg | BSA: m ² |
| Patient Comorbidities and Potential Incremental Risk Factors CABG Y <input type="checkbox"/> N <input type="checkbox"/> PPM Y <input type="checkbox"/> N <input type="checkbox"/> Porcelain Aorta Y <input type="checkbox"/> N <input type="checkbox"/> AAA Y <input type="checkbox"/> N <input type="checkbox"/> TAA Y <input type="checkbox"/> N <input type="checkbox"/> | | EOA needed to achieve an iEOA > 0.85 cm ² /m ² : LVEF %: Creatinine Clearance (cc/min): Mitral Regurgitation: Mild <input type="checkbox"/> Moderate <input type="checkbox"/> Severe <input type="checkbox"/> Previous BAV Y <input type="checkbox"/> N <input type="checkbox"/> Date: Pre-existing conduction disturbance: | | |
| Implant Date | TBD | (DD-MM-YYYY) | | |
| Comments | Patient appears to meet all sizing criteria for a 29mm Evolut valve. Bilateral transfemoral access appears ≥ 5.0mm. Suitable for 14F FX delivery system. Overlap = RAO: 13°, CAUD: 27° - Coplane = LAO: 10°, CAUD: 3° | | | |


MEDTRONIC ANALYSIS

ANNULUS

| | | | | | |
|-------------------------|-------|---|------------------|---|------|
| Diameter (mm) | 21.6 | x | 26.2 | , | 23.9 |
| | Min | | Max | | Mean |
| Perimeter (mm) | 75.1 | | , Derived Ø (mm) | | 23.9 |
| Area (mm ²) | 440.2 | | , Derived Ø (mm) | | 23.7 |

LVOT

| | | | | | |
|-------------------------|-------|---|------------------|---|------|
| Diameter (mm) | 18.7 | x | 25.7 | , | 22.2 |
| | Min | | Max | | Mean |
| Perimeter (mm) | 70.6 | | , Derived Ø (mm) | | 22.5 |
| Area (mm ²) | 375.9 | | , Derived Ø (mm) | | 21.9 |



| | | |
|------------------------------------|------|-----------|
| Max Ascending Aorta Diameter (mm) | 30.2 | |
| Sinotubular Junction Diameter (mm) | 27.3 | x 27.7 |
| | Min | Max |
| Sinus of Valsalva Diameter (mm) | 32.4 | 29.6 30.6 |
| | LCC | RCC NCC |
| Sinus of Valsalva Height (mm) | 22.6 | 20.9 21.8 |
| | LCC | RCC NCC |
| Coronary Ostia Height (mm) | 16.5 | 17.2 |
| | Left | Right |

RIGHT

| | | | |
|---------------------------|-----|---|-----|
| CIA Min Diameter (mm) | 9.0 | x | 9.6 |
| EIA Min Diameter (mm) | 7.9 | x | 8.3 |
| Femoral Min Diameter (mm) | 8.2 | x | 8.3 |

LEFT

| | | | |
|---------------------------|-----|---|-----|
| CIA Min Diameter (mm) | 9.2 | x | 9.6 |
| EIA Min Diameter (mm) | 7.5 | x | 7.7 |
| Femoral Min Diameter (mm) | 7.9 | x | 8.2 |


RIGHT

| | |
|------------------------------|---|
| Subclavian Min Diameter (mm) | x |
|------------------------------|---|

LEFT

| | |
|------------------------------|---|
| Subclavian Min Diameter (mm) | x |
|------------------------------|---|

Aortic Root Angle 47.9°



Please review images for direct aortic evaluation.

Calcium: Mild ☐ Moderate ☐ Severe ☐


VIV ADDITIONAL MEASUREMENTS

| | | |
|---------------------------------|--------|--------|
| Valve to Coronary Distance (mm) | To LCA | To RCA |
| Valve to STJ Distance (mm) | LCC | RCC |


Procedural Considerations

Aorta

ANNULUS




SOV DIAMETER





LVOT





STJ



Aorta




SINUS HEIGHT



AORTIC ROOT




IMPLANTER'S VIEW




Overlap

Aortic Root Angle



Coplane

Femoral Access - Right



RIA


Right Common Iliac Diameter
Min. Ø: 9.0 mm
Max. Ø: 9.6 mm
Avg. Ø: 9.3 mm

EIA (R)

Right External Iliac Diameter
Min. Ø: 7.9 mm
Max. Ø: 8.3 mm
Avg. Ø: 8.1 mm

RFA

Right Femoral Diameter
Min. Ø: 8.2 mm
Max. Ø: 8.3 mm
Avg. Ø: 8.2 mm




Femoral Access - Left



LIA

Left Common Iliac Diameter
Min. Ø: 9.2 mm
Max. Ø: 9.6 mm
Avg. Ø: 9.4 mm





EIA (L)

Left External Iliac Diameter
Min. Ø: 7.5 mm
Max. Ø: 7.7 mm
Avg. Ø: 7.6 mm







LFA

Left Femoral Diameter
Min. Ø: 7.9 mm
Max. Ø: 8.2 mm
Avg. Ø: 8.0 mm



Patient Evaluation Criteria

| Valve Size Selection | | Evolut™ PRO+ Bioprosthesis | | | |
|-----------------------------------|---------|---|---|---|---|
| | |  |  |  |  |
| Size | | 23 mm | 26 mm | 29 mm | 34 mm |
| Annulus Diameter | 23.9 mm | 17*/18-20 mm | 20-23 mm | 23-26 mm | 26-30 mm |
| Annulus Perimeter† | 75.1 mm | 53.4*/56.5-62.8 mm | 62.8-72.3 mm | 72.3-81.7 mm | 81.7-94.2 mm |
| Sinus of Valsalva Diameter (Mean) | 30.9 mm | ≥ 25 mm | ≥ 27 mm | ≥ 29 mm | ≥ 31 mm |
| Sinus of Valsalva Height (Mean) | 21.7 mm | ≥ 15 mm | ≥ 15 mm | ≥ 15 mm | ≥ 16 mm |
| Oversizing Percentage | | -4% | 9% | 21% | 42% |

*Measurement for TAV-in-SAV only. | †Annulus Perimeter = Annulus Diameter x π


NOTE: Evolut™ PRO+ valve size selection is identical to Evolut™ R valve size selection criteria

| Selection Criteria | |
|--|--|
| Access Consideration by MSCT | IFU Guidance by MSCT |
| Minimum Transarterial Access Vessel Diameter | Evolut PRO+ and Evolut R 23/26/29 TAVs ≥ 5.0 mm Evolut PRO+ 34 TAV ≥ 6.0 mm Evolut R 34 TAV ≥ 5.5 mm |
| Aortic Root Angulation. Femoral Access | Not recommended if >70 degrees. |
| Aortic Root Angulation. Left Subclavian | Not recommended if >70 degrees.** |
| Aortic Root Angulation. Right Subclavian | Not recommended if >30 degrees.** |
| Vascular Access Location. Direct Aortic Access | Ascending aorta access site ≥60 mm from basal plane.†† |

**Use caution in patients with a preexisting patent left internal mammary artery/right internal mammary artery (LIMA/RIMA) graft.

††For direct aortic access, ensure access site and trajectory are free of patent RIMA or preexisting patent RIMA graft.

Heart Team Procedure Plan (provided by hospital)

| | | | | | |
|--------------------------------------|--|---|--|--|--------------------------------|
| Patient Name | |  | | | |
| Planned Evolut™ R/ Evolut™ PRO+ Size | | <input type="checkbox"/> 23 mm | <input type="checkbox"/> 26 mm | <input type="checkbox"/> 29 mm | <input type="checkbox"/> 34 mm |
| Access Route | | <input type="checkbox"/> Transfemoral | <input type="checkbox"/> Subclavian | <input type="checkbox"/> Direct Aortic | |
| Vascular Access | | <input type="checkbox"/> Left <input type="checkbox"/> Right | <input type="checkbox"/> Left <input type="checkbox"/> Right | <input type="checkbox"/> Mini-Thoractomy | |
| | | <input type="checkbox"/> Percutaneous | | <input type="checkbox"/> Mini-Sternotomy | |
| | | <input type="checkbox"/> Cut-down | | | |
| Planned Anesthesia | | <input type="checkbox"/> General | <input type="checkbox"/> Conscious Sedation | | |

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. The Evolut™ R/ Evolut™ PRO+ transcatheter aortic valve has been approved by FDA for specific patient populations only.

Refer to the Instructions for Use for a full list of warnings, precautions, indications, and adverse events.

CAUTION: This report is provided pursuant to the terms of the Case Planner Physician Use Agreement and is based on information and images provided by the physician to Medtronic. This report is intended to be a resource to support physicians in their determination of proper case selection, device sizing and procedure planning, and is in no way intended to constitute medical advice or in any way replace the independent medical judgment of a trained and licensed physician with respect to any patient needs or circumstances. Physicians must conduct their own measurements and make their own medical judgments based on all of their patient's clinical and diagnostic records and images. Physician is solely responsible for all decisions and any medical judgments relating to patient diagnosis and treatment, including case selection and sizing of the device. Please see the complete Instructions of Use for all product indications, contraindications, precautions, warnings, and adverse events.

© 2013, 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo, and Further, Together are trademarks of Medtronic. ™* Third party brands are trademarks of their respective owners.

All other brands are trademarks of a Medtronic company.

Medtronic

710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Tel: (763) 514-4000
Fax: (763) 514-4879
Toll-free: (800) 328-2518

LifeLine
CardioVascular Technical Support
Tel: (877) 526-7890
Tel: (763) 526-7890
Fax: (763) 526-7888
rs.cstechsupport@medtronic.com

www.medtronic.com

Europe
Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale
CH-1131 Tolochenaz
Switzerland
Tel: (41 21) 802-7000

Canada
Medtronic of Canada Ltd
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
Tel: (905) 460-3800
Toll-free: (800) 268-5346

Asia
Medtronic International Ltd
49 Changi South Avenue 2
Nasaco Tech Centre
Singapore 486056
Singapore
Tel: (65) 6436-5000

Latin America
Doral Corporate Center II
3750 NW 87th Avenue, Suite 700
Miami, FL 33178
USA
Tel: (305) 500-9328

©2018 Medtronic.
Medtronic, Medtronic logo
and Further, Together are
trademarks of Medtronic.
All other brands are
trademarks
of a Medtronic company.
UC201402558a EN
07/2018